

## Technical Specifications

Nominal Size	Model FIXI	Model FOX
Weight (kg) / product	5.2	5.3 (without basket)
Max. body weight (kg)	100	100
Depth (mm)	625	670
Width (mm)	575	575
Adjustable handle height (mm)	600-690 / 5 positions	690-770 / 5 positions
Seat height (mm)	425	500
Seat width (mm)	385	385

## Spare Parts and Accessories

Art. No.	Name	Material
286.00.45	Wheel D=150 mm, black/yellow	Plastic
280.10.48	Padding for back bow	PUR
280.00.60	Brake set 1.265x1.000 Niro incl. nipple	VA steel
286.65.20	Seat for basket, black (for Fox)	
286.65.25	Seat without basket, black (for Fox)	
286.88.20	Seat for basket, green (for Fox)	
286.88.25	Seat without basket, green (for Fixi & Fox)	Pressed wood
189.10.10	Basket (only for "FOX")	Plastic encapsulated steel

Technical specifications subject to change without notice!



## User Instructions

BA No.: 02.05.012/07-07(KP)



**REBOTEC** Rehabilitationsmittel GmbH  
 D-49610 Quakenbrück, Artlandstr. 57-59  
 Telephone: +49-(0)5431/9416-0  
 Telefax: +49-(0)5431/9416-66  
<http://www.rebotec.de>  
 E-mail: info@rebotec.de

Serial No. / Product



## Rollator

Art. No. 186.00.10 „FIXI“



Art. No. 186.00.20 „FOX“



Also available with a basket

### Intended Use / User Advisory

The Rollators are designed for persons with limited walking ability and serve to increase their ability to move. The Rollators may be used only on solid, even surfaces.

The aluminium tubular design provides security against tipping over. Its lightweight not only makes the Rollator easier to drive, the bow handle also allows it to be pulled up towards the upper body, thus allowing the patient to easily overcome street kerbs. In addition the backrest may be outfitted with a back padding (accessory), which makes sitting comfortable. For transport, the Rollator can be folded up. The bow handle is also height adjustable to suit the height of the patient.

**Have your specialised dealer or care provider advise you on using the Rollator and its adjustment for proper function!**

## Instructions for Assembly and First-time Use

### 1. Unfolding the Rollator

Place the Rollator on its front wheels and push the frame with the rear wheels apart (A) **Fig. 1**. The Rollator clicks by itself into its final position (B) **Fig. 2**.

### 2. Adjusting the Handles

Next remove both star-grip screws (C) **Fig. 3**. Then choose the appropriate height (D) **Fig. 3** and adjust the grip handles to the next screw hole. Reinstall the star-grip screws (C) and tighten them firmly.

### 3. Attaching the basket (only for "FOX")

Using the 2 hooks, hang the basket (E) **Fig. 4** onto the tubular bow (F) of the Rollator.

### 4. Using the Brakes

By pulling up the left and right brake levers (G) **Fig. 5**, the back wheels are locked respectively. The permanent brake is set by pushing the brake levers (H) downwards. By pulling up the brake levers (G) the permanent brake is disengaged.

### 5. Folding the Rollator together

If present, first remove the basket (E). Then with one hand pull the **red cord** (J) **Fig. 2** upwards until the click-lock disengages. Then with the other hand lift the Rollator at the tubular bow (F).

### 6. Using the Seat (special accessory)

The seat (M) **Fig. 2** is screwed firmly to the frame (F).



- Before using the seat, activate the permanent brake (H), see Step 4.
- When the sun shines the seat may get hot! Therefore, check the temperature of the seat with your hand. If necessary, allow the seat to cool down, or place a cushion on the seat.

### 7. Adjusting the Brakes

Loosen the counter nuts (L) and then adjust the brake with the adjustment screw (K). Be sure that when the brake is released, the wheels are allowed to turn freely, and when the brake levers are pulled the wheels are firmly arrested. Then retighten the counter nut (L).



### Safety Instructions and Inspection

- After the first time set-up and before each use, always conduct a visual inspection of the Rollator for deficiencies. Check the screw connections to be sure they are tight. Tighten if necessary using the appropriate tools. If you have any doubt, then do not use the Rollator. Contact a specialised dealer.
- Before each use, check over the folding mechanism to be sure it has clicked into its final position.
- Check the function of the brakes and the permanent brake function as well.
- Use the grasp handle only for raising the Rollator and as a backrest.
- **DANGER OF TIPPING!**  
For moving the Rollator only use the rubberized handles.
- **DANGER OF TIPPING!**  
Never load the basket while using the Rollator with more than **4 kg** (9 lbs.)

## Care and Maintenance

The surface of the aluminium tubing has been provided with a paint coating. The Rollator can be cleaned with a commercially available cleaning or disinfecting agent – always follow instructions for use. All moving parts are permanently lubricated. **Inadequate or negligent care and maintenance leads to limitations or annulment of product liability by the Manufacturer.**

**Allow your Rollator to be checked over by specialised dealer once a year!**

### Advisory!

Always keep the User Instructions in case of future use.

Always bring a copy of the User Instructions with you whenever you leave the house so that in the case of an emergency another person will be able to operate the Rollator.

The identical serial number in the User Instructions is also printed on the nameplate of the product. For complaints or questions, always provide the model, serial no. and the User Instructions no.

### Warranty

We provide a warranty period of **24 months** for our Rollators. Wear parts are excluded.

### Re-Use (Rental)

Here we refer expressly to the strict conditions of the German product liability law (ProdHaftG §1 Abs.2 Nr.1).

### Figures:

